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EN DIRECTIONS FOR USE FOR WALL REINFORCEMENT IMPLANTS

CAREFULLY READ by the surgeon and the staff of the operating theatre and the staff concerned by the reception and the manipulation of implants. The use of this device has to be carried out exclusively in operating rooms with highly requested aseptic conditions.

The use of this device is strictly reserved for the specialists surgeons. The choice of implants must to be subject of a particular attention from the user by taking into account the encountered indication, the morphology, the activity, the physical and/or mental condition of the patient, this to limit the undesirable side-effects and the complications.

The results per and post-operative obtained by the user and the medical team cannot be correct in the case:

- The incorrect choice of various components of implants,
- A bad combination of various components of implants,
- A modification or a transformation of the device,
- An incorrect indication
- Defective technical operating
- A misunderstanding of the current standards, the commercial and/or technical documents (technical operating) (1), the device (2), the instruments,
- A care defect or a re-education defect after the intervention,
- Non compliance with the contents of this note,

And in these cases the user will be only held responsible (so loosening the responsibility of Microval).

The user of this device must be aware and must inform the patient before the operation,

- the precautions to be taken after the operation as well as the contra-indications and the undesirable side-effects which could arise
- that this device can allow a level of lesser activity,
- some activities which can be banned partially or totally in the time (3),
- that the inevitable wear of the device or / and the possible separation in the time can oblige the user or another surgeon to operate.

PRECAUTIONS OF USE : Before operation, please check that all specific instruments for the operation of devices are available and functional. Avoid any contact with objects which can damage the device. The damaged devices or / and having been in contact of a patient must be isolated and disinfected before cleaning and possible forwarding.

Mind, defect of fixation, positioning, assembling of the device can lead abnormal constraints and/or reduce the shelf life.

DESCRIPTION : These implants are made with knitted polypropylene, non-woven polypropylene or polyester. They could be enrolled and introduced through a trocar.

INDICATIONS : Abdominal wall reinforcement, Treatment of hernia and eventration.

CONTRA-INDICATIONS: Children during their growth. Intensive violent physical activities. Allergy reaction. Serious illness containing a risk of dangerous post-operative complications. The infection, the septicemia establish absolute contra-indications.

UNDESIRABLE SIED-EFFECTS, POTENTIAL COMLICATIONS:

POST-OPERATIVE : Wound haematoma. Slow wound healing. Tissue reactions as a result of allergy of the implanted material. Modification of the holding or the migration of the device. Infection. Adhesions. Seroma.

PRESERVATION : The device must be delivered and stored in its Microval original package, the storage conditions must avoid any risk of deterioration of the labelling, the package, and the device. The package must not be modified or transformed.

STERILITY : This device is delivered sterile and the sterilization method is mentioned on the outer package.

For any further information about the use of this device, please contact Microval.

- (1) Not detailed (for any further information please contact Microval)
- (2) If necessary, respect imperatively the side and/or the orientation of the device (see labelling)
- (3) Defined by the user

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