

## Disposable Balloon Trocars



### Features:

- + Universal seal, compatible for a wider size range of instruments.
- + Locking mechanism between obturator and cannula, provides a safe insertion.
- + Balloon and fixation design provides great stabilization. No damage to the abdominal wall.
- + Both air and saline can be input to the balloon.
- + Three kinds of tip available: Bladeless, Bladed and Optical.

Product Ref. Code	Product Name	Specification (mm)
CIT-BT-05	Disposable Balloon Trocar, Optical	5 x 100
CIT-BT-10		10 x 100
CIT-BT-12		12 x 100
CIT-BTW-05	Disposable Balloon Trocar, Bladeless	5 x 100
CIT-BTW-10		10 x 100
CIT-BTW-12		12 x 100
CIT-BTY-05	Disposable Balloon Trocar, Shield Bladed	5 x 100
CIT-BTY-10		10 x 100
CIT-BTY-12		12 x 100

CE 0197

# Disposable Hasson Trocar



## Features:

- + Durable universal seal, no reducer needed.
- + Recessed stopcock valve.
- + Unique fixing design.

Product Ref. Code	Product Name	Specification (mm)
CIT-HT-05	Disposable Hasson Trocar	5 x 100
CIT-HT-10		10 x 100
CIT-HT-12		12 x 100