


INTRODUCING **LAGIPOINT®**

For Effective Specimen Retrieval, Extracorporeal Anastomosis and Single Incision Surgeries!!

- Integrated seal system does not require trocar use, eliminating possible interference of trocar cannulas in the abdominal cavity
- Removable seal cap makes it easy to extract large specimens without removing the device
- Offers unmatched articulation, permitting use of standard laparoscopic instrumentation
- Maintains pneumoperitoneum for continuous access and visualization
- Accommodates 1.5cm to 3.5cm incisions and varied abdominal wall thickness for wide spectrum of procedures
- Reduces the risk of wound infections




Exclusively Distributed by:



Lap Tech Medical Sdn Bhd
Tel: +603 - 7880 7179
Website: www.laptechmedical.com.my
Email: admin@laptechmedical.com.my

Laparoscopic Cap



Removable Cap for Easy Specimen Retrieval

- Enables a laparoscopic approach before and after specimen retrieval
- Multi ports of 5mm, 12mm and 15mm seal enable use of a wide range of instrumentation
- Integrated seal system eliminates the need for trocars
- Accommodates straight, bent, angled, or curved instrumentation

360° Atraumatic Retraction & Wound Protection

- Distributes force evenly for optimal field of view
- Allows for maximum exposure with a minimum incision size
- Allows visualization of wound margins
- Soft materials and evenly dispersed tension promote patient post-operative comfort.

Greater Exposure

- 360° of elastic tension maximizes working area providing significantly greater exposure

Easily Placed and Removed

- Inserted in seconds. Flexible ring slips gently into incision

Significantly decreases wound infection

- Helps protect wound from contamination while helping to keep wound edges moist.

LAGIS LapShield™ Wound Protector/Retractor



Product Code	Number of Ports		Incision Length
	5mm	5-12mm	
TRL-0110R	1	1	1.5 – 3.5cm
TRL-0210R	2	1	
TRL-0220R	2	2	
TRL-0310R	3	1	
TRL80-0220R	2	2	3 - 6cm

Exclusively Distributed by:



Lap Tech Medical Sdn Bhd
Tel: +603 - 7880 7179
Website: www.laptechmedical.com.my
Email: admin@laptechmedical.com.my