

Exclusively Distributed by:



Lap Tech Medical Sdn Bhd  
Tel: +603-7880 7179  
Website: [www.laptechmedical.com.my](http://www.laptechmedical.com.my)  
Email: [admin@laptechmedical.com.my](mailto:admin@laptechmedical.com.my)



SIGMOID INSTRUMENT

## SIGMOID Suction Instrument

The Sigmoid clear, single-use Suction Instrument is primarily used for sigmoidoscopies, colonoscopies and femur suctioning during hip replacements to evacuate blood, fluid and other debris.

- Designed for use with a rigid sigmoidoscope
- ON/OFF suction control to reduce tissue grab
- Available with select tubing pre-attached
- 18 French, 30cm



POOLE INSTRUMENT

## POOLE Suction Instrument

Designed to evacuate pooled blood, fluid and debris during deep abdominal surgeries and cesarean sections, Poole single-use specialty suction instruments feature three versatile positions that quickly convert from gentle sump action to pinpoint suctioning capabilities.

- Removable sheath for pinpoint suctioning
- Reduce Tissue grab with removable sheath
- Removable sheath has 132 large eyes to provide efficient and rapid suction

### Ordering Information

Description	Case	Catalog Number
<b>Sigmoid Instrument</b>		
Sigmoid Suction Instrument	50	0033040
Sigmoid Suction Instrument w/6' (1.8m) Tubing pre-attached	20	0033050
<b>Poole Instrument</b>		
Poole Suction Instrument	50	0035040
Poole Suction Instrument w/6' (1.8m) Tubing pre-attached	20	0035050
Poole Suction Instrument w/10' (3m) Tubing pre-attached	20	0035070