

## FineRetriever<sup>™</sup>

Features a push/pull rod and an introducer tube to facilitate bag opening and tissue capture during laparoscopic procedures.

In the fully deployed condition, the bag opening is maintained in a fully-open position by a metallic rim.

A string with a slip knot enables the closing of the bag after the specimen collection.



- Strong, tear resistant
- Transparent pouch
- Radio-opaque rim



## SimpleBag

For quick, simple, safe specimen retrieval

PRODUCT CODE	DESCRIPTIONS
EBH-0507	5x7 320ml 10mm
EBH-0709	7x9 1000ml 12mm
EB-0507	5x7 320ml 10mm
EB-0709	7x9 1000ml 12mm

Package : 6 EA / Box