


- Non-slip handle
- Designed for easy penetration
- Available in standard and longer sizes
- Retractable safety, spring loaded tip protector
- Visual needle position indicator
- Clear handle for easy hang drop test
- High flow STOPCOCK



Lagis has designed insufflation needles that combine safety and high performance with an economical price.



Product Code	Description	Cannula's OD/ Length
IN-120	Disposable Insufflation Needle	2.1mm/120mm
IN-150	Disposable Insufflation Needle	2.1mm/150mm

*Packaging: 12/Box



Lap Tech Medical Sdn Bhd
Tel: +603-7880 7179

Website: www.laptechmedical.com.my
Email: admin@laptechmedical.com.my